EN

The **SISSSEL® Ball** with **SECUREMAX® system** offers you increased safety during use.

The SISSSEL® Exercise Ball is a medical device in accordance with EU Regulation No. 2017/745 and is suitable for use in sports, pilates, gyms and as a therapy in addition to the targeted balance and coordination training. It serves as an aid for strengthening the abdominal muscles, as well as the stabilizing muscles of the lower back. Furthermore, it enables the user to promote active and dynamic sitting at home or at work.

What is special about the SECUREMAX® system is the nature of the material used. This prevents the ball from bursting in the event of sudden, puncture damage (e.g. from sharp objects). In such cases, the ball loses its gradually, thereby allowing sufficient reaction time and reducing the risk of accidents. In order for the special properties of the SISSSEL® SECUREMAX® Exercise Ball to take effect and to ensure safety during use, the following must be observed. Please read the information carefully and keep the instructions for use for reference.

How to inflate the SISSSEL® SECUREMAX® Ball correctly:

1. Allow the flattened ball to warm to a good temperature for 24 hours before inflating. A room temperature of approx. 23 °C is ideal. Adequate warming makes inflation easier and prevents overinflation.

2. Inflate up to 80% of the specified ball diameter in several stages (please use the measuring tape provided). It is best to use a specially engineered hand pump (optional accessory) or a foot pump. 3. Leave to sit for at least 2 hours. After the specially designed pump machine has taken a little more time to adjust to the stretch and the internal pressure, it is ready to be pumped to the designated ball size (please use the measuring tape again for this) and firmness. 4. Never exceed the specified ball diameter. Due to the special material properties, the ball can initially feel a bit hard or still have slight lumps. This is normal and will go away after a short time. 5. After inflating, let the SISSSEL® SECUREMAX® Exercise Ball sit at room temperature for 24 hours before using it for the first time.

Please note: Once inflated, exercise balls gradually lose air. The diameter decreases by 2-3 cm per month. This is normal and is not a defect.

Important: The plug must be fully inserted into the exercise ball. Please use the attached plug remover to pull up the plug.

Safety instructions:

1. Do not use the ball if you suffer from pronounced sensory processing disorders or if you have acute lumbar spine syndrome, unstable fractures, acute lumbar spine injuries (e.g. herniated discs, inflammatory processes etc.), acute pain, paraesthesia or movement disorders.

1. Where possible, use the ball under the instruction and supervision of trained specialists. 2. Avoid overly difficult exercises or high-risk experimentation. 3. Do not use the ball on uneven or pointed edges (e.g. keys, scissors, pens, etc.). 4. Regularly check the integrity of the ball surface. 5. Do not use the ball if it is damaged. 6. Keep the ball away from radiators or other heat sources and avoid direct sunlight. 7. The ball should not come into direct contact with painted or coated surfaces. 8. Use the ball for a dry surface only.

General information:

1. The shape of the ball may change slightly when used. After being left for a while, the ball will return to its original shape.

1. In order to clean the ball, wipe with a damp cloth. Then rub dry. 2. Distinct using regular surface cleaning. 3. The ball must only be filled with air. 4. Dispose of damaged balls correctly and never repair them. 5. The SISSSEL® SECUREMAX® Ball is designed for a maximum static weight limit of 500 kg. The maximum user weight is 130 kg.

1. Burst resistant up to 130 kg in the event of puncture damage caused by a pointed object. 2. The SISSSEL® SECUREMAX® Exercise Ball is made of phthalate-free PVC and is recyclable.

Important information, please note:
The SISSSEL® SECUREMAX® Exercise Ball is made with CIPTROFOL® BI. CIPTROFOL® is an excellent alternative for products with high demands. These are products for which standard plasticizers (phthalates) should not be used either on the grounds of banned substances or in accordance with the precautionary principle. Your new ball may therefore smell a vinegar or vinegar cleaner (instead of a chemical smell), but this will go away completely after a few days.

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