

FOR USE WITH



INTELECT RPW 2
Falcon Handpiece



MOBILE 2 RPW
D-Actor Handpiece



RPW LITE
Sparrow Handpiece

C15



- Ø 15 mm, CERAmA-x® transmitter
- Penetration depth: 0 – 30 mm
- Any type of enthesopathy

R040 BEAM



- Ø 15 mm, AkuSt transmitter
- Penetration depth: 0 – 35 mm
- RPW Lite version is black
- Any type of enthesopathy

D15 GOLDEN DEPTH



- Ø 15 mm, trigger point transmitter
- Penetration depth: 0 – 50 mm
- Deep-lying pathologies, chronic conditions

D20-S STANDARD OSCILLATOR



- Ø 20 mm, AkuSt transmitter
- Penetration depth: 0 – 50 mm
- Muscle tissue and connective tissue

D20-T GOLDEN OSCILLATOR



- Ø 20 mm, trigger point transmitter
- Penetration depth: 0 – 50 mm
- Muscle tissue and connective tissue: high-energy application

FOR USE WITH



INTELECT RPW 2
Falcon Handpiece



MOBILE 2 RPW
D-Actor Handpiece

A6 POINT



- Ø 6 mm, AkuSt transmitter
- Penetration depth: 0 – 20 mm
- Shock wave acupuncture

T10 FINGER



- Ø 10 mm, trigger point transmitter
- Penetration depth: 0 – 40 mm
- Trigger points, deep-lying enthesopathy

F15 WHITE SOFT



- Ø 15 mm, CERAmA-x® transmitter
- Penetration depth: 0 – 30 mm
- Superficial pain regions, muscles of mastication

B15



- Ø 15 mm, CERAmA-x® transmitter
- Penetration depth: 0 – 30 mm
- Superficial pain regions, sensitive areas

R15



- Ø 15 mm, AkuSt transmitter
- Penetration depth: 0 – 35 mm
- RPW Lite version is black
- Any type of enthesopathy

D35-S LARGE OSCILLATOR



- Ø 35 mm, CERAmA-x® transmitter
- Penetration depth: 0 – 60 mm
- Muscle tissue and connective tissue, large muscles

V10



- Ø 10mm
- Treatment of specific muscle points
- Compatible with V-Actor II and V-Actor HF

V25



- Ø 25mm
- Muscle smoothing
- Massage
- Increase in local blood circulation
- Compatible with V-Actor II and V-Actor HF

V40



- Ø 40mm
- For large Muscle areas
- Muscle smoothing
- Massage
- Increase in local blood circulation
- Compatible with V-Actor II and V-Actor HF

FOR USE WITH



INTELECT RPW 2
Falcon Handpiece

ATLAS



- Soft flexible tipped transmitter 15mm for highly sensitive regions
- Any type of enthesopathy
- Penetration depth 0-15mm

SPINE-ACTORS



SPAN: 20 mm

SPAN: 27.5 mm

SPAN: 35 mm

- Double Tip shape to allows treatments of the paravertebral region.
- They come in three sizes allowing treatment in all sections of the spinal column muscles (lumbar, thoracic and cervical).

FASCIA TREATMENTS



PERI-ACTOR KNUCKLE
Treatment of all medium-sized and large fascia



PERI-ACTOR SCOOP
Treatment at or under bone edges



PERI-ACTOR SCRAPER
Treatment of all fascia of the musculoskeletal system



PERI-ACTOR SPHERE
Spot treatment of fascia

- Unique shape especially designed to dissolve clots in the Fascia structures

RPW TRANSMITTERS